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CARE		STAGES	ACCESSO- RIES	MODES	TIME / PRESSURE LEVEL	CARE		STAGES	ACCESSO- RIES	MODES	TIME / PRES- SURE LEVEL	
SCARS	Retractable Adhesive Fibrous	1. Mobilize	Suction head of your choice	Pulsed	5min / Level 2 (150 mBar)	SPINAL COLUMN	Neck pain	1. Vascularize	R15 / R25 suction head	Conti- nuous	5min / Level 3 (200 mBar)	
		2. Smooth		Conti- nuous				2. Mobilize	Suction head of your choice	Pulsed		
	Caesarean sections	1. Vascularize	R35 suction head	Conti- nuous	5min / Level 1 (100 mBar)			3. Soften		Conti- nuous		
		2. Smooth	Palpate and roll		10min / Level 3 (200 mBar)			Back pain	1. Vascularize	Suction head of your choice		Conti- nuous
DRAINAGE	Edema	1. Mobilize	Suction head of your choice	Pulsed	5min / Level 3 (200 mBar)		2. Mobilize		Pulsed			
		2. Drain		Conti- nuous	10min / Level 3 (200 mBar)		3. Soften		Conti- nuous			
	Upper limbs	1. Stimulate	Suction head of your choice	Pulsed	5min / Level 3 (200 mBar)		Low back pain		1. Vascularize		Suction head of your choice	Conti- nuous
		2. Mobilize						2. Mobilize	Pulsed			
		3. Drain		Conti- nuous				3. Soften	Conti- nuous			
	Lower limbs	1. Stimulate	Suction head of your choice	Pulsed	5min / Level 3 (200 mBar)		Sciatica	1. Vascularize	Suction head of your choice	Conti- nuous	5min / Level 3 (200 mBar)	
		2. Mobilize						2. Mobilize		Pulsed		
		3. Drain		Conti- nuous				3. Soften		Conti- nuous		
TENDINOPATHY		1. Relieve	Suction head of your choice	Pulsed	5min / Level 3 (200 mBar)	BURNS		1. Mobilize	Suction head / Palpate and roll	Pulsed	5min / Level 2 (150 mBar)	
		2. Mobilize		Conti- nuous				2. Soften	Suction head of your choice	Conti- nuous	5min / Level 3 (200 mBar)	